


Welcome to Teledent.com.au Dental Imaging Diagnostics

How the site works:

STEP ONE: log in

From Teledent's home page (www.teledent.com.au), you can access your secure Patient and Sessions Work page

Enter your email address and password.



Step two: check progress of your reports or upload a new

case Once logged in, you will have a choice of:

- seeing the progress status of previous cases submitted previously to Teledent
- download the prepared radiology reports
- resend receipts for your financial records
- submitting a new case for reporting
- changing your password



Teledent

Specialist Dentomaxillofacial
Radiology Reporting

GET A REPORT

Show My Reports

Get a new report

Welcome back Louise

You are signed in as
louise@perio.net.au

[Change User?](#)

Under Show my Reports:

CHANGE MY PASSWORD

Reports in progress


	Louise Brown cc Panoramic report	2015-12-09 - Order:000010	In Progress
	jojo abd Panoramic report	2015-12-10 - Order:000013	In Progress
	testing Accutoma data Cone beam report Medium and Small FOV	2015-12-19 - Order:000018	In Progress

Completed reports ready for downloading

	Louise Brown Test case 1 Cone beam report large field of view FOV	2015-12-08 - Order:000009 Total: \$ 0.00 GST: \$ 0.00	Download: CBCT report Oral Diagnostics Clinic 59189 PDF.pdf	Re-Send Invoice
---	---	---	--	-----------------

Step three:

Submit a new case for reporting



Click on Get a new report and the choice of reports will show:



Teledent

Specialist Dentomaxillofacial
Radiology Reporting

Our Reports




Cone beam report

\$ 200.00

Comprehensive radiology report will include interrogation and analysis of entire dataset for pathology; a reconstructed panoramic view and one or two other views to show relevant findings.

[Download sample report](#)

[GET REPORT](#)



Intraoral film

\$ 70.00

A radiology opinion and report on intraoral images (bitewing radiograph, periapical radiograph and/or occlusal images)

[Download sample report](#)

[GET REPORT](#)



Lateral Cephalograms Report

\$ 65.00

A radiology report on a standard 2D lateral cephalograms image

[GET REPORT](#)



Panoramic report

\$ 90.00

A report on a standard 2D panoramic image

[Download sample report](#)

[GET REPORT](#)

Select the report type (this example shows small/medium FOV cone beam) and

- enter your name
- enter the required patient information
- it is important to complete the clinical information, including identifying the primary reason for requesting the report. The radiologist will always address the region of interest and clinical purpose firstly, and then complete a full pathology review.
- zip the dicom datafile - the folder **MUST** be zipped so that it can be opened at the other end (see separate instructions)

You can upload the files directly to Teledent by clicking on the Upload button in the Patient Files box (please be aware the larger the file the longer it takes to upload the file. Size of file and your internet speed determine the amount of time).

Acceptable file types are multi-file dicom (dcm) files, or Invivo (inv) files. Multi-file DICOM folders **MUST** be compressed (zipped) prior to uploading.

For OPGs and cephs- please upload a HIGH resolution jpeg image or the dicom data

- click on box to browse to find patient datafile or drag file into the top box

- you can add additional images, previous reports, clinical photos into the lower box. Any digital file can be included with a case such as Word documents, JPEGs,

Report Information

Referring Physician Name

Patient Name

D.O.B Gender Age

Patient Address

Clinical Notes (Please be comprehensive)

Patient Image/Dataset

Drop file here or Click to Upload

Relevant Previous Images/Datasets

Drop file here or Click to Upload



Cone beam report
Comprehensive radiology report will include interrogation and analysis of entire dataset for pathology; a reconstructed panoramic view and one or two other views to show relevant findings.

\$ 200.00

All prices exclude GST



SECURE PAYMENTS
We accept the following cards:





- it is essential that the CBCT folder has been fully uploaded before you proceed to the next step

Once upload is complete, scroll down this same page to either

Additional Services

- Work up of portfolio of images from CBCT dataset \$ 150.00

Include multiple slices for implant site measurements, impacted teeth including identification of IAN and full work up of TMJ/airways/paranasal sinuses cases. PDF of 1:1 sized worked up images can be requested

- RUSH Report \$ 150.00

For any reporting required with a 24 hour turnaround

 Add To Cart

I Agree to Terms of Use

[Terms of use](#)

MAKE PAYMENT >

Keep me in the loop



- Once you have entered your credit card details, the final step that you need to do is to click the orange PLACE ORDER NOW button



- You will receive an email confirming the order, and a receipt for the payment.
- Once the report has been completed, you will then receive an email asking you to log onto your secure site to download the report



Under Show my Reports:

The FAQ section has a lot more information – you will find most of your queries will be answered if you look through this section of the website:

- A: How does the Teledent site work?
- B: How do I export the cone beam data files to upload into Teledent?
- C: What is a dentomaxillofacial radiologist?
- D: Some medical radiologists specialise in dental radiology?
- E: What is the usual turnaround time for image reports?
- F: Is there a Medicare or private health insurance rebate for Teledent radiology reporting?
- G: Do you provide surgical guides for implant cases?
- H: Can patients send Teledent their scans or speak to a radiologist?